

#### PRECISION IS OUR PROFESSION





## Agenda

MeKo MedTech

Vasculoy® - Ni-free Alloy for Permanent Stents

Questions & Answers

## Where the story begins







In May 1991, Clemens Meyer-Kobbe started MeKo as a one-man business.





MeKo Manufacturing e.K.


Laser material processing and post processing

More than 300 qualified employees

Development of new materials, optimizing of material properties

Dedicated to high precision and challenging processes

High export rate with > 50 % outside Europe





## Huge Stock for Rapid Prototyping



- More than 3.000 different tube and sheet lots on stock
- Available materials:
  - > 316L medical, L605, Phynox, MP35N
  - > NiTi (Nitinol)
  - > Bioresorbables: Mg / Resoloy<sup>®</sup>, Fe, Zn, polymers / PolymediX<sup>®</sup>
  - > Vasculoy® (Nickel-free alloy)
- Securing material quality thanks to in-house inspection and measuring instruments for material properties







## Manufacturing Processes

for Medical Devices





#### Laser Cutting / Drilling

- Nearly no limitations in cutting length and tube diameter
- Superior dimensional accuracy

#### Laser Welding

- Spot welding and seam welding
- In-process quality control

#### Heat Treatment / Shape Setting

- Extremely fine grain sizes
- Controlled heating cycles ensure reliable Af-temperature setting

#### Electropolishing

 Widely adjustable range of material removal

#### Passivation and Final Cleaning

- Automated passivation / cleaning line for consistent results
- Enhanced biocompatibility

#### Quality Inspection

- 100 % inspection of all implants
- Visual and automatic microscopes





## PRECISION IS OUR PROFESSION

\_MeKo.de

Vasculoy®

The new and only Ni-free alloy for permanent stents





## Agenda



- 1) Why a Ni-free Alloy?
- 2 Material Requirements for Stents
- 3 Vasculoy® The New Material?
- 4 Material Properties of Vasculoy®
- 5 Vasculoy® in comparison
- 6 Road Map of a Vasculoy® Stent



## 1 Allergy



- Studies show that approximately **13** % of the population are allergic to nickel with a strong upward trend [1].
- Individuals, who are allergic to nickel, demonstrate a greater restenosis rate [2, 3].
- Furthermore, the corrosion products of 316L are toxic to the primary culture of vascular smooth muscle **cells** when the nickel concentration is higher than 11.7 ppm [4].





# 1 Allergy



Patients with a <u>severe</u> Ni-allergy are seldom but show immediately significant symptoms.

Two samples are shown in the next slides.





# 1 Allergy



Pruritic and an eczematous dermatitis worse on the right leg after NiTi stent implantation [5]





Resolution of dermatitis one month following stent explantation [5]



# Allergy





Location of iliac stent (arrow) [6]



A picture on admission showing rash around abdomen [6]

#### After iliac stent explantation:

"Postoperatively, her dermatitis significantly improved, and her depression disappeared immediately." [6]



# Regulatory Requirements



- The European Parliament and Council passed the Nickel Directive:
  - > limits **nickel content** in products used for epithelization after piercing to 0,5 % mass of Ni to total mass [7]
  - > limits **nickel release** from objects intended for use in direct and prolonged contact with skin to 0,5 μg/cm² per week [7]
- The FDA released two information regarding implants with Ni-content:
  - Biological Responses to Metal Implants [8]
  - Guidance for Industry and Food and Drug Administration Staff: "Corrosion of implant devices made of, or having components made of, nickel-rich alloys (e.g., nitinol, stainless steel, MP35N) results in the release of nickel ions, which may lead to various modes of toxicities." [9]




# Responsibility and Liability



Who takes the **responsibility** for implanting a nickel-containing device if an identical Ni-free device is available?

Can Ni-allergic patients claim for **indemnification** if they experience restenosis?

Patients must be **tested for Ni-allergy** before implantation of a Ni-containing implant.









#### Requirements Material properties Mechanical properties Strength Thin strut sizes Strain / Plasticity Highly deformable without cracks Corrosion resistance No degradation or ion releases Biocompatibility No negative tissue reactions, Ni-freedom MRI compatibility No artifacts under MRI Radiopacity Angiographic visibility under X-ray

## 3 Vasculoy® - The New Material?




- There are many Ni-free <u>industrial</u> alloys:
   all have a high strength but **no** plasticity → not suitable for stents
- A special alloy for vascular applications was developed by MeKo:
  - Based on stainless steel but Nickel is replaced by Nitrogen
    - > nitrogen stabilizes the austenite grain structure
    - > with the nitrogen content in the alloy the material properties are adjusted for stents




= Vascular Alloy



# 4 Strength and Strain







*Vasculoy*®

- High strength and break elongation
- Nearly identical to L605




## 4 Recoil



Tests performed according to ISO 25539-2 and ASTM F 2079

	316L			Vasculoy®			CoCr - L605		
Stent Location	Distal	Middle	Proximal	Distal	Middle	Proximal	Distal	Middle	Proximal
Diameter <sub>inflated</sub> in mm	3,87	3,87	3,88	3,76	3,8	3,87	3,79	3,79	3,83
Diameter <sub>final</sub> in mm	3,72	3,72	3,75	3,56	3,6	3,63	3,58	3,58	3,6
Stent Recoil in %	4,0	3,5	3,5	5,5	5,5	6,0	5,5	5,5	6,0





Recoil equivalent to L605




### (4)

## **Corrosion Resistance**



Cyclic potentiodynamic analysis according to ASTM F2129 is used to study the corrosion behavior of metallic materials



	316L	Vasculoy®	CoCr - L605
<b>Protection Potential</b> E <sub>p</sub> in mV	859	1040	591
<b>Breakdown Potential</b> E <sub>b</sub> in mV	823	842	769
Passive range in mV	825	1058	802



Best corrosion resistance of all stent materials



# Biocompatibility



Biocompatibility investigations according to EN ISO 10993

	Vasculoy®
Cytotoxicity test (EN ISO 10993-5)	<b>No</b> proliferation inhibition
Hemolysis elution method (EN ISO 10993-4)	<b>No</b> hemolytic effects
Hemocompatibility dynamic conditions (EN ISO 10993-1)	<b>No</b> activation of the coagulation system
Chemical analysis quantification of inorganic leachable substances (EN ISO 10993-1)	<b>No</b> quantifiable leachable substances
Chemical analysis detection of organic substances (EN ISO 10993-12)	<b>No</b> detectable substances



- Perfect biocompatibilityNo compromises




### (4)

# **MRI** Compatibility



Tests were performed according to the standard test method ASTM F2052, ASTM F2119 and ASTM F2182

	316L	Vasculoy®	CoCr – L605
Magnetically Displacement Force 3 Tesla Weight of Probe in Nx10 <sup>-4</sup> Displacement Force in Nx10 <sup>-4</sup> Force resulting deflection angle in °	2,65 0,6 13	2,88 0,3 7	2,79 0,24 5
Radio Frequency Induced Heating 1,5 Tesla Probe / Reference Probe in °C	0,8 / 0,9	0,7 / 0,8	1,3 / 1,1
MR Image Artifacts 1,5 Tesla Spin / Gradient echo artifact in mm	4,7 / 6,4	2,6 / 5.5	3,4 / 4,6
MR Image Artifacts 3,0 Tesla Spin / Gradient echo artifact in mm	5,5 / 8,2	4,4 / 6,4	3,8 / 6,0



MRI 3.0 T gradient echo image with transversal alignment of stent: L605 (top), Vasculoy (bottom).

*Vasculoy*°


- Results superior to 316L
- Results comparable to CoCr L605



# 4 Radiopacity

\_\_MEDTECH

- Vasculoy® has a similar density/radiopacity like 316L
  - > Recommendation for radiopaque markers
    - Two on each stent side.
    - No compromises to stent material
    - Good visibility of stent markers (no artifacts)
- Radiopaque materials like PtCr show following disadvantages:
  - High costs
  - Lower strength and strain
  - Less corrosion resistance
  - Artifacts in high resolution angiography (no insight view into stent!)



Angiographical visibility of implanted 316L stent without markers (left) and of a similar stent with markers (right).



## Comparison



Mechanical Properties:

•	Strength	316L	<	Vasculoy®	=	L605
---	----------	------	---	-----------	---	------

- Strain / Plasticity 316L < Vasculoy® = L605
- Corrosion resistance 316L < Vasculoy® >> L605
- Biocompatibility
   Vasculoy® is the only Ni-free stent material
- MRI compatibility 316L < Vasculoy® = L605
- Radiopacity
   Radiopaque markers enable good visibility

*Vasculoy*°

- The best material for stents and other implants
- The only choice for nickel allergic patients




# 6 Road Map of a Vasculoy® Stent

# \_\_\_\_MEDTECH

#### MeKo's customer support for the development

- Customer sends a design identical to an existing L605 stent
- MeKo manufactures samples / prototypes (cost < 2000 EUR)</li>
- Initial testing of Vasculoy® stents by customer
  - > results will be equivalent to the existing L605 stent
- Customer decides for Vasculoy® stent registration:
  - > all necessary biocompatibility data concerning the material will be provided by MeKo for free
  - > documentation support with lots of Ni-allergy literature
- Serial production at MeKo for prices like L605 stents





# Vasculoy® List of Literature



- [1] Thomas, P., Schuh, A., Ring, J. et al. Orthopädisch-chirurgische Implantate und Allergien. Orthopäde 37, 75–88 (2008). doi: 10.1007/s00132-007-1183-3.
- [2] Köster, R., Vieluf, D., Kiehn, M. et al. Nickel and molybdenum contact allergies in patients with coronary in-stent restenosis. Lancet. 2000 Dec 2; 356(9245):1895-7. doi: 10.1016/S0140-6736(00)03262-1.
- [3] Schmidt M., Goebeler M. Nickel allergies: paying the Toll for innate immunity. J Mol Med (Berl). 2011 Oct; 89(10):961-70. doi: 10.1007/s00109-011-0780-0.
- [4] Shih, C.C., Shih, C.M., Chen, Y.L. et al. Growth inhibition of cultured smooth muscle cells by corrosion products of 316 L stainless steel wire. J Biomed Mater Res. 2001 Nov;57(2):200-7. doi: 10.1002/1097-4636(200111)57:2<200::aid-jbm1159>3.0.co;2-5.
- [5] Jetty, P., Jayaram, S., Veinot, J. et al. Superficial femoral artery nitinol stent in a patient with nickel allergy. J Vasc Surg. 2013 Nov;58(5):1388-90. doi: 10.1016/j.jvs.2013.01.041.
- [6] Guntani, A., Kawakubo, E., Yoshiga, R. et al. Metallic allergy requiring removal of iliac stent: report of a case. Surg Case Rep. 2020 Apr 26;6(1):82. doi: 10.1186/s40792-020-00843-8.
- [7] European Union. European Parliament and Council Directive 94/27/EC. Official Journal of the European Communities. 1994, June 30<sup>th</sup>. https://eur-lex.europa.eu/legal-content/EN/TXT/PDF/?uri=CELEX:31994L0027
- [8] U.S. Food & Drug Administration. Biological Responses to Metal Implants. Sept. 2019. https://www.fda.gov/media/131150/download?attachment
- [9] U.S. Food & Drug Administration. Select Updates for Non-Clinical Engineering Tests and Recommended Labeling for Intravascular Stents and Associated Delivery Systems. 2015, Aug. 15<sup>th</sup>. https://www.fda.gov/media/93134/download







Two stents are on the table, the design and mechanical properties are identical, one with Nickel, the other without:



Which stent will you take?



